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I. CONTENTS

1 syringe + 2 protective tubes for intraligament anaesthetic syringes.

II. DESCRIPTION OF SYMBOLS USED

DANGER SIGNALS

WARNING :
important
information !

Please wear
goggles

Please wear gloves

REMARKS

General indications
with no danger to
people or objects

Visual
examination

Please read this !

III. AREA OF APPLICATION

The syringe is used for dental requirements in the area of the field of local or regional anaesthesia, any other form of use is forbidden and can cause danger.

 For this syringe, we applied Community Directive 93/42/CEE, transposed by the French regulations on medical equipment (Official Bulletin - L169 - 12th of July, 1993).

 In accordance with these provisions, the contra angle must be used by someone with experience of dental medicine, for the application described, and in compliance with the provisions in effect in relation to the prevention of accidents at work, protection of employment and the instructions for use.

Preparation and maintenance of this medical appliance have to be done exclusively by someone having training on prevention against infections, on self-protection and on protection of the patients.

In accordance with these provisions, the user must:

- only use working cartridges and needles without defects,
- observe the principle of correct use,
- protect themselves and the patient or third party against any danger,
- prevent all contamination by the product.

IV. GENERAL SAFETY INFORMATION

Before use, check that the syringe is not damaged and that no part of it is missing.

No modification or addition may be made to the device without the company's express prior agreement.

Significant risk of injury and contamination when handling the device and needles. We recommend using a needle recapper (Anthogyr, ref. 2470).

Remove the needle immediately after the anaesthetic and place it in a suitable collector. Comply with current legislation on disposal of medical waste.

Risk of complications related to the anaesthetic: The user must have been trained in first aid and have suitable resuscitation equipment available. When choosing an anaesthetic, the patient's specific counter-indications must be taken into account (carry out an appropriate medical questionnaire).

Choosing the type of needle, the type of anaesthetic, the quantity injected and the speed of the injection are the user's responsibility. The user is also responsible for the aspiration technique (particularly in the case of block or para-apical anaesthetics).

All anaesthetic cartridges which have been used, even partially, must be thrown away.

Risk of the anaesthetic cartridge bursting: use a protective tube, particularly for intraligament or intraseptal anaesthetics.

V. TECHNICAL SPECIFICATIONS

The cartridge capacity, the type of thread on the needle and, as appropriate, the type of aspiration, are specified on the product label.

VI. USING THE SYRINGE

Risks owing to non-sterile products:
The syringe is not sterile when supplied.

The syringe must be sterilised before it is first used. (See paragraph on "Hygiene and maintenance").

Before use, check that the syringe is not damaged and that there are no parts missing.


Only use high quality anaesthetic cartridges which conform to ISO11499.

Use a syringe which is suitable for the anaesthetic cartridge's capacity. See label.

Only use high quality needles which conform to ISO7889.

For needles with a non-metric thread on the hub, use a syringe with a reference including the letters US.

6.1. Syringe with hinge :


1. Pull back the rod of the piston thoroughly and hold it.

2. Rock the tube of the syringe to open the back hole of the tube.

3. Insert the injection cartridge, the sealed part first.

4. Put the tube of the syringe in its initial position holding the rod pulled.

6.2. Syringe with lateral insertion of the cartridge :

1. Pull back the rod of the piston thoroughly.

2. Insert laterally the injection cartridge, the sealed part first, in the tube of the syringe.

6.3. Syringe with manual aspiration system (type 2a) :

1. Turn the handle clockwise sticking the corkscrew into the piston plug of the cartridge or push the handle to insert the harpoon into the piston plug of the injection cartridge.

2. Aspiration will be applied with a back movement of the piston.

6.4. Self aspirating syringe (type 2b) :

Release the injection pressure to create aspiration within the cartridge.

It may happen that the self-aspirating principle does not satisfactory work with some cartridges or some needles. You must then test the right work of the self-aspiration before the first use.

6.5. Bayonet syringe (intraligament):

1. Unscrew the bayonet and remove the syringe barrel.
2. Insert the cartridge, capsule part first, into the back part of the syringe barrel.
3. Screw on the syringe barrel by engaging the bayonet.

Screw the needle following the instructions given by needles manufacturer.¹

¹: Plug with standard metric screwing except for references followed by US.

VII. HYGIENE AND MAINTENANCE **in accordance with standard ISO 17664**

Re-sterilization of reused medical appliances has to be done by someone who has been trained and is protected, and the regulations in force have to be adhered to. The re-sterilization protocol must be adapted to the infectious risk.

For each used product : please refer to the manufacturer's instructions. Respect in particular the concentrations, the exposure durations, the solutions' replacement and the products' lifetime. Never mix the products. Please dispose of products as per their instructions.

Wear adapted protective clothes.

In order to avoid any risk of infection and injury, it is vital to wear protective gloves.

Plastic protective tubes for intraligament syringes are not suitable for autoclave sterilisation: they must be replaced after every use.

7.1. Pre-disinfection

Pre-disinfection has to be made immediately after the operation. The anaesthetic cartridge and the needle must be removed from the syringe.

• By soaking :

Use only detergent or disinfectant solutions with a neutral or slightly alkaline pH. It is extremely inadvisable to use products which may fix X proteins to the unit (eg. alcohol, aldehydes...).

Carefully rinse and dry.

• Or : Thermo disinfection thermal disinfection for equipment meeting the requirements of standard ISO 15883, immediately after the intervention:

The thermo-disinfection cycle lasts 10 minutes at 93°C (203°F). Follow the manufacturer's instructions for use.

Use only low alkaline cleaning products pH < 9.

Ensure the drying cycle is completed.

Dry any water residue with compressed air :
Wear glasses and a mask when you do this.

7.2. Cleaning

Use only detergent or disinfectant solutions with a neutral or slightly alkaline pH. It is extremely inadvisable to use products which may fix X proteins to the unit (eg. alcohol, aldehydes...).

Do not use sodium hypochlorite : There is a high risk of corrosion.

• **Brushing / Immersion :**

Brush meticulously every part and surface with a soft brush.

Allow the items to soak for the correct period of time.

Carefully rinse and dry each part.

> **Or : Cleaning in an ultrasonic Bath :**

Only use products which are compatible with this technique. Keep the parts separate from each other as contact in the bath will cause unsightly external marks.

7.3. Operating test

This test has to be done prior to each sterilization.

Check that the piston rod slides freely.

7.4. Sterilization

The syringe **must be sterilized before the first use and after each use.**

The syringe can be sterilized at 135°C (275°F) **only** in a vacuum steam autoclave for 20 minutes minimum (sterilization maintenance period).

We strongly advise to use class B autoclaves.

Any other sterilization method must not be used.

Read the autoclave manufacturer's instructions. Respect the space between each bag and do not overload the autoclave.

- Before being sterilized, the instruments must be pre-disinfected, cleaned, lubricated and tested.
- Remove the anaesthetic and the needle before sterilisation.
- Sterilization pouches must be adapted to the syringe and to the autoclave and in accordance with the standard NF EN868.

Only one single syringe must be put in a pouch.

Keep the syringe and its accessories within the sterilization bag in a clean and dry place.

Warning : After each sterilization, check that there is no residual water.

VIII. REPAIR

In the event of breakdown, please contact your approved distributor or the Anthogyr after sales service department directly.

All repairs must be carried out with parts and subassemblies approved by the manufacturer.

Repairs must only be carried out by an approved distributor or by the after sales service department at the factory.

For any revision or repair the instrument must be returned complete and sterile with proof of sterility. It must be sent together with a document stating the problem as well as the complete name and address of the dentist.

For any claims to be considered under the warranty, please enclose a copy of the invoice or the delivery note with the equipment.

Replacement of parts is guaranteed for 7 years after the product ceases to be marketed.

AFTER-SALES DEPARTMENT CONTACT DETAILS

Services S.A.V.
ANTHOGRYR
2237, Avenue André Lasquin
74700 Sallanches - FRANCE
Tel direct: +33 (0)4 50 58 50 53
Mail: sav@anthogyr.com

IX. ACCESSORIES

To be ordered from your approved distributor:

<i>Description</i>	<i>Ref.</i>
Aspiration corkscrew	3804 F
Batch of 12 protective plastic tubes for 1.8cc intraligament syringes	2429
Batch of 12 protective plastic tubes for 2.2cc intraligament syringes	2429/22
Needle recapper	2470

X. WARRANTY

This equipment is guaranteed for parts and labor against any defect in manufacture for a period of 12 months from the date of invoicing. This warranty does not apply to wearing parts and does not cover transport costs.

For claims to be considered under the guarantee, please enclose a copy of the invoice or the delivery note with the equipment.

Any modification or addition to the product without Anthogyr's express agreement will nullify this guarantee.

The guarantee is invalid in the event that the technical instructions supplied with the equipment are not followed.

Anthogyr cannot be held responsible for damage and its consequences resulting or possibly resulting from normal wear and tear, incorrect use, cleaning or maintenance, non-compliance with the conditions for use or connection, scaling or corrosion, impurities in the air or water supply systems, or abnormal chemical or electrical effects or non-compliance with the Anthogyr instructions for use, maintenance or installation and any other instructions by the manufacturer.

XII. DISPOSAL

Respect current legislation.

In particular: used needles must be placed in a suitable collector.

Conform to current legislation regarding the disposal of medical waste.

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