


# Anthofit® HE angulated conical abutments

The Anthofit® HE Regular system (R platform) is increasing the list of prosthetic components with the introduction of new angulated conical abutments of 18° and 30°.









→ Ref. HE PA001  
Price : please contact your local consultant  
Angulated conical abutment 18° - Height 4 mm



→ Ref. HE PA002  
Price : please contact your local consultant  
Angulated conical abutment 30° - Height 4.5 mm

**These angulated conical abutments allow angle correction required for screw-retained restorations (bar-overdenture prosthesis or bridges). They add to the existing range of straight conical abutments and allow exhaustive use of this type of prosthetic restorations.**

These 2 new references are compatible with the secondary parts below (already mentioned in the catalog):

HE IMPLANT	R
<b>Healing screw</b>	The angulated conical abutments are not supplied with a healing screw. During the healing stage, place the protective cap on the abutment.
<b>Pop-in transfer</b>	 HETR008
<b>Protective cap</b>	 OIC0022
<b>Analog</b>	 HEAN006
<b>Temporary abutment</b>	 OIC0023
<b>Castable cylinder</b>	 OIC0024
<b>Gold cylinder</b>	 OIC0025