

Rotating angulated conical abutments


oxiom[®] angulated conical abutments are now available with **360°** rotation!



Increased flexibility for multi-unit implant restorations!

These abutments allow for maximum flexibility of positioning in the patient's mouth (360° indexing possibilities).

The impression must be taken at the abutment level in order to maintain perfect accuracy for the restoration. Impression taken at the implant level does not allow for accurate repositioning of the abutment in the mouth.



→ 360° rotating angulated conical abutments
Price: please contact your local consultant

→ Ref. OPACR22 - 360° rotating angulated conical abutment - Ht. 2,5 mm - ang. 18°

→ Ref. OPACR23 - 360° rotating angulated conical abutment - Ht. 2,5 mm - ang. 30°

→ Ref. OPACR32 - 360° rotating angulated conical abutment - Ht. 3,5 mm - ang. 18°

→ Ref. OPACR33 - 360° rotating angulated conical abutment - Ht. 3,5 mm - ang. 30°