

## LABSIDE compatibility list

Our libraries are compatible with CAD software : Exocad®, 3Shape® et DentalWings®.

They can be downloaded from [www.anthogyr.com](http://www.anthogyr.com) section DIGITAL / LABSIDE.

Before any installation, it is important to identify :

- . The CAD software and its version,
- . The component references for Axiom® BL and Axiom® TL
- . The scanning accessories

Before any installation, it is important to uninstall previous versions of libraries that are no longer useful.

V - 26/06/2020

Implant Platform		Milling options		Workflow for conventionnel impression		Workflow for digital impression with IOS (Intra-Oral Scanner)		Exocad® All versions			3Shape® Scanners and CAD software			DentalWings Scanners and CAD software
Components references	Components type	Indexed	Rotational	Scan-Adapters for lab scanners	Analogues for plaster model	Digital transfers for IOS	Analogues for printed model	Lab scanner Imetric®	Lab scanner Medit®	Others configurations				
								version 5 and upper	version 2.0.0.0 and upper	all opened Exocad® versions	version 2014 and upper	versions 2017 and 2018	version 2019 and upper	version 4.0.0 and upper
Name of the associated library files														
OPFLEX413/423/433 OPFLEX513/523/533	FlexiBase® Axiom BL	✓	-	152-27-SAA 152-27-SAA	OPIA100 OPIA100	152-27-DT 152-27-DT	152-27-PA 152-27-PA	Anthogyr_Labside_imetric_V2.02.01	Anthogyr_Labside_exoCAD_V2.02.01			Anthogyr_Labsidemb_2017_3shape_V2.02.01	Anthogyr_Labsidemb_2019_3shape_V2.02.01	Anthogyr_Labsidemb_DentalWings_V2.02.01
MUNFLEX MUFLEX	FlexiBase® Multi-Unit	-	✓	151-04-SAO 151-03-SAO	MUNA100 MUA100	151-04-DT-MUN 151-03-DT-MU	151-04-PA 151-03-PA							
TFLEX-N TFLEX-R	FlexiBase® Axiom TL	✓	-	156-0X-SAA 156-0X-SAA	TA100-N TA100-R	156-01-DT 156-02-DT	156-01-PA 156-02-PA							
TFLEX-N-P TFLEX-R-P		-	✓	156-0X-SAO 156-0X-SAO	TA100-N TA100-R	156-01-DT-IL 156-02-DT-IL	156-01-PA 156-02-PA							
OPFLEX403 OPFLEX503	FlexiBase® Axiom BL	✓	-	152-27-SAA 152-27-SAA	OPIA100 OPIA100	-	-	exocad_imetric		Anthogyr flexibase 3shape	-	-	dentalwings	